

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-7751

Bay-Ker, Up, RR, Cer90°, 180x5mm



- Ceramic
- Original
- Offset Up
- Shaft Length 180mm
- Chip Set Ejector
- 40° Angled Thin Foot Plate

Bay-Ker, Up, RR, Cer90°, 180x5mm

[http://buxtonbio.com/baykeruprrcer90180x5mm\\_1680](http://buxtonbio.com/baykeruprrcer90180x5mm_1680)