

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-8734

Bay-Ker, L-hnd, RR, Cer180x3, 40°, tfp



- Ceramic
- Original
- Offset Left Handed
- Shaft Length 180mm
- Chip Set Ejector
- 40° Angled Thin Foot Plate

Bay-Ker, L-hnd, RR, Cer180x3, 40°, tfp

[http://buxtonbio.com/baykerlhndrrcer180x340tfp\\_1709](http://buxtonbio.com/baykerlhndrrcer180x340tfp_1709)