

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

86-4051

BAP Frs-Sm Punch, 180x5mm, 90°up



- TAP Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 180 mm Shaft Length

TAP Ferris-Smith Punch, 180x5mm,
90° up-angled

http://buxtonbio.com/bap_frsm_punch180x5mm90up_1164