

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 86-4462

BAP Frs-Sm Punch, 230x6mm, 40° up



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 230 mm Shaft Length

TAP Ferris-Smith Punch, 230x6mm,  
40° up-angled

[http://buxtonbio.com/bap\\_frsm\\_punch230x6mm40up\\_1122](http://buxtonbio.com/bap_frsm_punch230x6mm40up_1122)