

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

863-4962

BAP F-S Pnch, 280x6, 40° up, Ejec



Photo  
Pending

- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Ferris Smith Handle
- 280 mm Shaft Length

TAP Frs-Sm Punch, 280x6mm, 40°  
up-angle, Ejec

[http://buxtonbio.com/bap\\_fs\\_pnch280x640upejec\\_1140](http://buxtonbio.com/bap_fs_pnch280x640upejec_1140)