

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

863-7622

TAP, Bay-Kers'n, Up, Ej, 40°, 150x2



Take-Apart Small Kerrison Handle
Punch Rongeur, Up Bayonet, 150mm
length, 2.0 bite size, 40° foot angle, Up
biting, standard foot plate

http://buxtonbio.com/tapbaykersnupej40150x2_1569