

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 863-8154

TAP, Bay-Ker, RH, Ej, 110x5, 40° tfp



Take-Apart Small Kerrison Handle  
Punch Rongeur, Right-Handled  
Bayonet, 110mm length, 5.0 bite size,  
40° foot angle, Up biting, thin foot plate

[http://buxtonbio.com/tapbaykerrhej110x540tfp\\_1619](http://buxtonbio.com/tapbaykerrhej110x540tfp_1619)