

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 863-8352

TAP, Bay-Ker, R-hnd, Ej, 180x5, 40°



Take-Apart Small Kerrison Handle  
Punch Rongeur, Right-Handled  
Bayonet, 180mm length, 5.0 bite size,  
40° foot angle, Up biting, standard foot  
plate

[http://buxtonbio.com/tapbaykerrhndej180x540\\_1634](http://buxtonbio.com/tapbaykerrhndej180x540_1634)