

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

863-8534

TAP, Bay-Ker, LH, Ej, 110x3, 40°fp



Take-Apart Small Kerrison Handle  
Punch Rongeur, Left-Handled  
Bayonet, 110mm length, 3.0 bite size,  
40° foot angle, Up biting, thin foot plate

[http://buxtonbio.com/tapbaykerlhej110x340fp\\_1589](http://buxtonbio.com/tapbaykerlhej110x340fp_1589)