

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-4251

BAP Fr-Sm Pnch, Cer, 200x5, 90° up



- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 200 mm Shaft Length

TAP Frs-Sm Punch, Cer, 200x5mm,  
90° up-angle

[http://buxtonbio.com/bap\\_frsm\\_pnchcer200x590up\\_1399](http://buxtonbio.com/bap_frsm_pnchcer200x590up_1399)