

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

867-5222

BAPSmH Pnch, Cer, 20x2, 40° up, TFP



Photo  
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 200 mm Shaft Length

TAP Large Handle Punch, Cer,  
200x2mm, 40° up, TFP

[http://buxtonbio.com/bapsmh\\_pnchcer20x240uptfp\\_1381](http://buxtonbio.com/bapsmh_pnchcer20x240uptfp_1381)