

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-5442

BAP Pnch, Cer, 230x4mm, 40° up, TFP



Photo
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 230 mm Shaft Length

TAP Sm-Hndl Punch, Cer,
230x4mm, 40° up, TFP

http://buxtonbio.com/bap_pnchcer230x4mm40uptfp_1382