

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-7531

HD, Cer Bay-Kers'nUp, 110x3, 90°



- Ceramic
- TAP
- Up Offset
- Shaft Length 110mm
- Chip Set None
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, HD Take-A-Part, Ceramic coated, 110 mm shaft, offset up, 90°, 3mm bite

[http://buxtonbio.com/hdcer\\_baykersnup110x390\\_1747](http://buxtonbio.com/hdcer_baykersnup110x390_1747)