

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-8214

BAP Pnch, CerR-H-B, 150x1, 40°TFP



Photo
Pending

- Ceramic
- TAP
- Right Handed Offset
- Shaft Length 150mm
- Chip Set None
- 90° Angled Foot Plate

TAP Punch, CerRt-Handl Baynt,
150x1mm, 40°TFP

http://buxtonbio.com/bap_pnchcerrhb150x140tfp_1813