

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-8614

HD Bay-Ker:L,Cer,40°tfp,150x1



- Ceramic
- TAP
- Left Handed Offset
- Shaft Length 150mm
- Chip Set None
- 40° Angled Thin Foot Plate

Bayonet Kerrison rongeur, HD Take-A-Part, Ceramic coated, 150 mm shaft, left-handed, 45° , 1mm bite, thin footplate

http://buxtonbio.com/hd_baykerlcer40tfp150x1_1785