

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-8744

BAP Pnch, Cer, LHB, 180x4, 40°TFP



- Ceramic
- TAP
- Left Handed Offset
- Shaft Length 180mm
- Chip Set None
- 40° Angled Thin Foot Plate

TAP Punch, Cer, Lf-Handl Baynt,  
180x4mm, 40°TFP

[http://buxtonbio.com/bap\\_pnchcerlhb180x440tfp\\_1793](http://buxtonbio.com/bap_pnchcerlhb180x440tfp_1793)