

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-3451

BAP Pnch, Cer, 230x5, 90° up, Ejec



Photo
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Small Handle
- 230 mm Shaft Size

TAP Sm-Hndl Punch, Cer,
230x5mm, 90° up, Ejec

http://buxtonbio.com/bap_pnchcer230x590upejec_1412