

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 868-3462

BAP Pnch, Cer, 230x6, 40° up, Ejec



Photo  
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Lumbar
- Chip Handling Ejector
- Small Handle
- 230 mm Shaft Length

TAP Sm-Hndl Punch, Cer,  
230x6mm, 40° up, Ejec

[http://buxtonbio.com/bap\\_pnchcer230x640upejec\\_1366](http://buxtonbio.com/bap_pnchcer230x640upejec_1366)