

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-5452

BAPPnch,Cer,230x5,40°up,Ej,TFP



- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling Ejector
- 230 mm Shaft Length

TAP Sm-Hdl Punch,Cer, 230x5mm,40°
up,Ejec, TFP

http://buxtonbio.com/bappnchcer230x540upejtfp_1385