

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-76-5102

TAP Pnch, Cer, U-Bay, 150x2, LCUp

Photo
Pending

- Ceramic
- TAP
- 180 mm Shaft Length
- Regular Foot Plate

Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur, Up
Bayonet, 180mm length, 2.0 bite size,
40° foot angle, Up biting, standard foot
plate

http://buxtonbio.com/tap_pnchcerubay150x2lcup_1848