

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-7664

TAP, BayKer, UpEj, Cer40° 150x6tfp



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur, Up
Bayonet, 150mm length, 6.0 bite size,
40° foot angle, Up biting, standard foot
plate

http://buxtonbio.com/tapbaykerupejcer40150x6tfp_1760