

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone:973-560-4848
Fax: 973-560-4999

868-77-5102

BAP Pnch,Cer,U-Bay,180x2,LCUp



Photo
Pending

- Ceramic
- TAP
- 180 mm Shaft Length
- Regular Foot Plate

Take-Apart Small Kerrison Handle ,
Ceramic CoatedPunch Rongeur, Up
Bayonet, 180mm length, 2.0 bite size,
40° foot angle, Up biting, standard foot
plate

http://buxtonbio.com/bap_pnchcerubay180x2lcup_1850