

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-77-5163

BAPnch,Cer,U-Bay,18x3,LCUp,TFP



Photo
Pending

- Ceramic
- TAP
- 180 mm Shaft Length
- Thin Foot Plate

Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur, Up
Bayonet, 180mm length, 3.0 bite size,
40° foot angle, Up biting, thin foot plate

http://buxtonbio.com/bapnchcerubay18x3lcuptfp_1851