

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-7751

TAP, Bay-Ker, Up, Ej, Cer90°, 180x5



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur, Up
Bayonet, 180mm length, 5.0 bite size,
90° foot angle, Up biting, standard foot
plate

http://buxtonbio.com/tapbaykerupejcer90180x5_1766