

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 868-8132

TAP, Bay-Ker, RH, Ej, Cer110x3, 40°



Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur,  
Right-Handled Bayonet, 110mm  
length, 3.0 bite size, 40° foot angle, Up  
biting, standard foot plate

[http://buxtonbio.com/tapbaykerrhejcer110x340\\_1807](http://buxtonbio.com/tapbaykerrhejcer110x340_1807)