

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-8324

TAP, Bayonet, Right offset, Ejector, Ceramic Coat, 180x2, 40°



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur,
Right-Handled Bayonet, 180mm
length, 2.0 bite size, 40° foot angle, Up
biting, thin foot plate

http://buxtonbio.com/tap_bayonet_right_offset_ejector_ceramic_coat180x2_40_thin_footplate_1825