

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-8341

TAP, Bayonet, Right offset, Ejector, Ceramic Coat, 180x4, 90°



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur,
Right-Handled Bayonet, 180mm
length, 4.0 bite size, 90° foot angle, Up
biting, standard foot plate

http://buxtonbio.com/tap_bayonet_right_offset_ejector_ceramic_coat180x4_90_1823