

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-8562

TAP, Bay-Ker, LH, Ej, Cer 110x6, 40°



Photo  
Pending

Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Left-  
Handled Bayonet, 110mm length, 6.0  
bite size, 40° foot angle, Up biting,  
standard foot plate

[http://buxtonbio.com/tapbaykerljejcer110x640\\_1779](http://buxtonbio.com/tapbaykerljejcer110x640_1779)