

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 868-8764

TAP, BaKer, LH, Ej, Cr180x6, 40° tfp



Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Left-  
Handled Bayonet, 180mm length, 6.0  
bite size, 40° foot angle, Up biting, thin  
foot plate

[http://buxtonbio.com/tapbakerlhejcr180x640tfp\\_1796](http://buxtonbio.com/tapbakerlhejcr180x640tfp_1796)