

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

190-121

Fiberoptic Rod, 150 x 6mm



Illuminate the wound sites with this flexible, fiber optic rod. An anodized aluminum fixture accommodates this rod and mounts on the rack of the Buxton Shoulder retractor or most any device with parallel surfaces up to 8 mm thick.

http://buxtonbio.com/fiberoptic_rod_150_x_6mm_10407