




## 24-6015

Mod.Kölbel Retractor, dbl prng



- For primary retraction
- Blunt tips
- Double Pronged
- Part of Standard and Deluxe Shoulder Sets.
- Part of our new Optimal Shoulder Set.

Used most often as primary retraction following initial incisions. Ring-handled, self-retaining retractor with long sweeping arms. This version features angled arms and curved pairs of blunted prongs. Prongs are arranged in parallel series at the distal portions of the delicately tapered arms and point laterally in the direction of the retraction force.