



## 24-6101

### Kolbel Retractor Set w/4-blade



- Set features our most popular frame, the Kölbel retractor
- Includes one pair of 36x36 mm blades
- Includes one pair of 36x53 mm blades
- Components of Standard and Deluxe Shoulder sets (24-6001 + 24-8001)

One-handed placement, positioning and operation make the Kölbel Retractor Sets well suited for a broad range of open and mini-open shoulder repair procedures. Spring loaded ball bearings set into each arm of the retractor frame hold blades while allowing limitless free rotation.