



## 24-6102

### Kölbel Retractor Frame only



- Our most popular Retractor extends 75 mm (with blades attached)
- Original frame; features a slight bend in the arms
- Modular; fits any of our Shoulder blades
- Available bundled with a variety of blades
- 200 mm overall frame length
- Components of Standard and Deluxe Shoulder sets (24-6001 + 24-8001)

One-handed placement, positioning and operation make the Kölbel frame well suited for a broad range of open and mini-open shoulder repair procedures. Spring loaded ball bearings set into each arm of the retractor frame hold blades while allowing limitless free rotation.