



24-6162

Modified Kolbel Retr frame



- Modified Kölbl retractor frame features a wingnut as a center screw
- Attach one of three middle arms to the frame using the wingnut
- Snap on any of our wide variety of side blades to the end hubs
- Remove the center blade if not needed during surgery
- Available in sets, bundled with different blade assortments
- 200 mm overall frame length
- Component of Optimal Shoulder set (24-8500)

The Kölbl frame is modified to accommodate swivel, snap-in blades at the outer reaches of its arms and slotted center blades designed to slide along the post extension to the pivot screw around which the arms articulate. See also items 24-6160, 24-6161, 24-6168 and 24-6169 for various frame and blade combinations you may wish to consider.