



24-6168

Modified Kölbel Retractor Set, with 3 middle blades + 6 Kölbel



- Modified Kölbel retractor frame features a wingnut as a center screw
- Attach one of three included middle arms to the frame using the wingnut
- Middle blades depths are 36 mm, 53 mm and 68 mm.
- Includes one pair of 36x36 mm blades
- Includes one pair of 36x53 mm blades
- Includes one pair of 36x68 mm blades
- All items are components of Optimal Shoulder set (24-8500)

The Kölbel frame is modified to accommodate swivel, snap-in blades at the outer reaches of its arms and slotted center blades designed to slide along the post extension to the pivot screw around which the arms articulate. The 24-6168 set includes the modified frame, 3 pairs of swivel blades 36 mm wide, in 3 lengths (36, 53 and 68 mm), and 3 sliding center blades to match each of the swivel blade lengths. See also 24-6162 and 24-6160 for frame and frame with third arm combinations.