



24-6179

Middle-Arm Upgrade w/3 Blades



- Upgrade is done at our NJ facility with quick same-day turn around. Get retractor back in as little as 48 hours!
- Contact us first by phone or email to initiate the upgrade.
- Clean certificate and RMA # is required before sending the upgrade in.
- Clean certificate is available in the "Records" dept of our website.

Convert your existing Kölbel retractor frame to accommodate an interchangeable center blade that will provide better exposure and stability. Upgrade includes the wingnut center screw conversion as well as three 20 mm wide interchangeable center arms, matching the depths of the popular 36 mm, 53 mm and 68 mm long swivel blades.