

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone:973-560-4848  
Fax: 973-560-4999

## 52-5257

Sq Hndl Hook,shrp,4-prong,7mm



- Overall length: 175 mm
- Shaft length: 75 mm
- 7 mm reach
- 4-Prong
- Sharp

Square Handle Hook, sharp, 4-prong,  
7mm

[http://buxtonbio.com/sq\\_hndl\\_hookshrp4prong7mm\\_612](http://buxtonbio.com/sq_hndl_hookshrp4prong7mm_612)