



## 52-5853

### Swivel Blade, 42 x 15 mm, with blunt teeth for Distal Radius



- Quick, ball snap connection to frame
- Swivel connection automatically distributes tissue tension
- Broad width enhances access and exposure
- Blunt prongs prevent sharp injuries
- Blade removing forceps (74-3199) available



Distal Radius Fracture Retractor  
Blade, 42x15 mm, with blunt teeth

<http://buxtonbio.com>

[com/swivel\\_blade\\_42\\_x\\_15\\_mm\\_with\\_blunt\\_teeth\\_for\\_distal\\_radius\\_fractures\\_retractor](http://buxtonbio.com/swivel_blade_42_x_15_mm_with_blunt_teeth_for_distal_radius_fractures_retractor)