



74-3117

Caspr R&P Retr.w/hinge&10bldes



- Caspar Rack & Pinion Retractor Set
- Hinged Rack & Pinion Frame, parallel excursion
- Frame maximum opening: 70 mm
- 5 Smooth Blades + 5 Short Tooth Blades
- Smooth Blade depths in mm: 45, 50, 55, 60, 65
- Toothed Blade depths in mm: 40, 45, 50, 55, 60
- Blade width: 20 mm



Caspar Rack & Pinion Retractor
w/hinged arms & 10 blades