

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane 8  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 74-3435

### Cloward Lamina Spreader



- Classic Cloward
- Length 180 mm
- Diverging jaws 4 mm wide
- Spreads to 17 mm

Cloward Lamina Spreader

[http://buxtonbio.com/cloward\\_lamina\\_spreader\\_537](http://buxtonbio.com/cloward_lamina_spreader_537)