

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 827-5301

IVD Rngr, Cer, up 30°, 180x1.5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Up 30° Cup Angle
- 180 mm Shaft Length

IVD Rongeur, Ceramic coated, up 30°,  
180x1.5mm

[http://buxtonbio.com/ivd\\_rngrcerup\\_30\\_180x1\\_5mm\\_812](http://buxtonbio.com/ivd_rngrcerup_30_180x1_5mm_812)