

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 83-6522

Caspar IVD Rongeur, dn30°, 140x2mm



- Standard IVD
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 140 mm Shaft Length

Caspar IVD Rongeur, down 30°,  
140x2mm, fenestrated cups

[http://buxtonbio.com/casparivdrongeurdn30140x2mm\\_872](http://buxtonbio.com/casparivdrongeurdn30140x2mm_872)