

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 83-7655

Caspar IVD Rongeur, straight, 300x5mm, fenestrated cups



- Standard IVD
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Straight Cup Angle
- 300 mm Shaft Length

Caspar IVD Rongeur, straight,  
300x5mm, fenestrated cups