

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 83-7683

Caspar IVD Rongeur, down angled, 300x3mm, fenestrated cups



- Standard IVD
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 300 mm Shaft Length

Caspar IVD Rongeur, down angled,  
300x3mm, fenestrated cups

[http://buxtonbio.com/caspar\\_ivd\\_rongeur\\_down\\_angled\\_300x3mm\\_fenestrated\\_cups\\_2003](http://buxtonbio.com/caspar_ivd_rongeur_down_angled_300x3mm_fenestrated_cups_2003)