

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 83-7686

Caspar IVD Rongeur, dwn, 300x6mm



- Standard IVD
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 300 mm Shaft Length

Caspar IVD Rongeur, down angled,  
300x6mm, fenestrated cups

[http://buxtonbio.com/caspar\\_ivd\\_rongeurdwn300x6mm\\_2003](http://buxtonbio.com/caspar_ivd_rongeurdwn300x6mm_2003)