

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

831-112-15

IVD Rngr, R-hndl, 110x1.5,up30



- Specialty IVD
- Specialty Bayonet
- Offset Right
- Up 30° Cup Angle
- Shaft Length 110

IVD Rongeur, right handled, 110mm
shaft, 1.5mm cup, up 30°

http://buxtonbio.com/ivd_rngr_rhndl_110x1_5up30_942