

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 831-114-30

IVD Rngr, R-hndl, 110x3mm, dn30



- Specialty IVD
- Specialty Bayonet
- Offset Right
- Down 30° Cup Angle
- Shaft Length 110

IVD Rongeur, right handled, 110mm  
shaft, 3mm cup, down 30°

[http://buxtonbio.com/ivd\\_rngr\\_rhndl\\_110x3mmdn30\\_945](http://buxtonbio.com/ivd_rngr_rhndl_110x3mmdn30_945)