

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 831-134-30

IVD Rngr, R-hndl, 130x3mm, dn30



- Specialty IVD
- Specialty Bayonet
- Offset Right
- Down 30° Cup Angle
- Shaft Length 130

IVD Rongeur, right handled, 130mm  
shaft, 3mm cup, down 30°

[http://buxtonbio.com/ivd\\_rngr\\_rhndl\\_130x3mmdn30\\_946](http://buxtonbio.com/ivd_rngr_rhndl_130x3mmdn30_946)