

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-6552

Caspar IVD Rngr, Cer, dn30° 160x2mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 160 mm Shaft Length

Caspar IVD Rongeur, Ceramic coated,
down 30°, 160x2mm, fenestrated cups

http://buxtonbio.com/casparivdrngrcerdn30160x2mm_900