

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-7552

Weil-Bsly Rngr, Cer, str, 300x4mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Straight Cup Angle
- 300 mm Shaft Length

Weil-Blakesley Rongeurs, Ceramic coated, straight, 300 x 4 mm

http://buxtonbio.com/weilbsly_rngrcerstr300x4mm_887